PIC/S GMP対応
Compatible with PIC/S GMP

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- Introducing the HTST heat treatment system used for the production process of active pharmaceutical ingredients for biopharmaceuticals as well as application for treating injectable solutions, eye drop solutions, etc.
This system is suitable for the inactivation of viruses in biopharmaceutical ingredients as well as the bioburden management of injectable solutions and eye drop solutions.

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Advantage

Capable of handling small lot and small flow rate

The small-diameter (6mm) double-tube type heat exchanger enables small capacity, HTST heat treatment.

Continuous treatment of approximately 20 liters per batch is possible.

High-yield

High yield has been achieved by use of the liquid substitution system equipped with a conical tank and an electric conductivity meter.

Compact design

Small footprint capable of fitting into available space has been made possible by the small flow rate specifications (50L/H) and ultra-compact specifications (1600Wx900Dx2200H).

Utilities


An electrical heat source can be selected for the small flow rate type device, and the device can be installed in an environment with no boiler facility.

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

Application examples

- ・HTST heat treatment system for biopharmaceutical ingredients (medium solution, etc.), injectable solutions, eye drop solutions, etc.

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Specifications

Performance	50L/Hr ~ 8000L/Hr	
Temperature specifications	25°C ~ 110°C (HTU) ~ 25°C	Temperature and holding time can be set according to items.
Target item	Low-viscosity fluid (medium solutions, injectable solutions, eye drop solutions, etc.)	
Heat exchange system	Indirect heating and cooling system by double tube sheet type shell and tube, or double tube	
Heat retention method	Holding tube (vacuum insulation)	
Heat medium	Hot-water heating system (electrical or steam) or steam heating system	
Piping bore diameter	6mm・12mm・15A・1S・1.5S・2S	Optimal flow rate is ensured according to each bore diameter.
Device/Sensor	Diaphragm valve, sensor, high-sanitation sensor	Compatible with bifurcation point 4D
	Rotary pump	
Material of wetted part	SUS316L-equivalent electrolytic grinding (+ passivation treatment)	

※ After discussing use conditions, we will determine the optimal tank capacity, heat exchanger, sensors, bore diameter, and material. (Testing using our company's test equipment is also available.)

※ This system is compatible with an automatic CIP system according to alkali and acid.

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IWAI PHARMA TECH CO.,LTD.